GENOSS Inflator

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A Product description

The GENOSS Inflator is an inflation/deflation device with an analog pressure gauge that registers pressures from VAC to 30 atm ospheres (VAC to 441 psi).

The GENOSS Inflator is consisted of a threaded plunger and locking switch which allow for the generation and sustaining of pressure. A three way stopcock should be included for use during preparation of the device and in conjunction with interventional devices.

B Indication for use

The GENOSS Inflator is recommended for use during vascular procedures in conjunction with interventional devices such as balloon dilation catheters to create and monitor pressure in the balloon.

C Warning

Use only with radiopaque inflation medium recommended by the dilation catheter manufacturer. The in vivo dilation catheter inflation pressure should not exceed the maximum inflation pressure recommended by the manufacturer of the dilation catheter. This device is intended for single use only; Do not reuse. Do not resterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.

D Direction for use

- 1) Connect the GENOSS Inflator to a 10 or 20cc syringe filled with contrast medium or mixed contrast medium with normal saline.
- 2) To fill the GENOSS Inflator with contrast medium:
 - a. Open the stopcock and rotate the locking switch counter clockwise to its UNLOCK position. Ensure that located on the locking switch is centered within the syringe portion of the GENOSS Inflator.
 - b. Aspirate contrast medium by slowly withdrawing the plunger.
 - c. Evacuate air in the GENOSS Inflator by slowly injecting the necessary volume from the GENOSS Inflator back into the syringe.
 - d. Repeat steps 2) b through 2) c until all the air is removed from the GENOSS Inflator and about 6cc of contrast remains.
- 3) To create a closed system, rotate the plunger locking switch clockwise into its LOCK position.
- 4) Connect the GENOSS Inflator to the inflation port of the interventional device such as balloon catheter hub.
- 5) When balloon inflation is desired, verify that the plunger locking switch is in its LOCK position. Rotate the plunger handle clockwise to increase pressure or counterclockwise to decrease pressure.
- 6) Pressure can be quickly released by rotating the locking switch one half turn counterclockwise, to the UNLOCK position and then retracting the plunger handle.
To maintain this position, turn the locking switch one half turn clockwise to the LOCK position.

- 7) To quickly release from the vacuum position, turn the locking switch one half turn counterclockwise to the UNLOCK position.

E Precautions

Note the "Use By" date specified on the package.

Inspect all products prior to use. Do not use if the package is open or damaged, or if product is damaged.

Refer to the appropriate device labeling for indications, contraindications, and potential complications associated with the use of interventional and/or diagnostic devices.

This device should be used by physicians trained in vascular procedures that utilize interventional and/or diagnostic devices such as balloon dilation catheters.

F Handling Cautions


This device is sterilized using Ethylene Oxide gas and not pyrogenic material. And also, this device is designed and intended for single use only. Do NOT reuse, reprocess or resterilize. GENOSS will not be responsible for any direct, incidental or consequential damages resulting from re-sterilization and/or reuse. Do not use if sterile barrier is damaged. If damage is found, call your GENOSS representative. After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy.

Do not expose the GENOSS Inflator to organic solvents, e.g., alcohol. If strong resistance is met during manipulation, stop the procedure and determine the cause of the resistance before proceeding.

G Storage

Store at room temperature (1°C ~30°C) in dry and dark places.

H Symbols

	Catalogue number		Batch code
	Do not reuse		Date of manufacture
	Caution		Use by
	Sterilized using ethylene oxide		Manufacturer
	Do not resterilize		Keep away from sunlight
	Fragile, handle with care		Do not use if package is damaged
	Consult instructions for use		Temperature limitation
	Quantity per package		

MS-P451/IFU-DS1102(Rev.15,2402)