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Hyaluronic Acid Dermal Filler with Lidocaine

MONALISA *Lidocaine*

Device Generic Name Sterile Absorbable Injectable Dermal Fillers

Description MONALISA Lidocaine is a sterile, transparent gel of stabilized hyaluronic acid of non-animal origin with 0.3% lidocaine. MONALISA Lidocaine is supplied in a glass syringe with a Luer-lock fitting. The contents of the syringe have been sterilized using moist heat and pressure. The product is for single use only. Disposable sterile needles are provided with syringe. Information about size of the needle is printed on its packaging.

Intended Use

MONALISA Lidocaine is intended to be used for soft tissue augmentation and restoration of facial lipotrophy and contain the local anesthetic lidocaine for a more comfortable injection in patients over the age 21.

Indication For Use

1. Precaution before Use

- Use the product after examining the packaging for damage or contamination of the product and checking the expiration date. Do not use expired, damaged or contaminated products.
- The patient must be informed of the purpose, expected results, precautions and possible adverse events before treatment.
- Proper injection technique is essential to obtain favorable results. Injection technique of the manufacturer prior to initial injection is important for a more successful injection.
- Only qualified medical professionals must administer the treatment.

2. Treatment Procedures

- The use of sterile needles and cannula matching the Luer-lock specifications of the syringe is essential for safe handling.
- The needle must be correctly attached to the syringe for safe handling. Incorrect attachment can result in the detachment of the syringe and the needle. Hold the junction of the syringe and the Luer-lock with one hand and the sterile needle shield with the other. Rotate the needle to securely fix it to the syringe.
- The treatment area is to be disinfected with aseptic solution.
- Do not bend the needle to prevent breakage.
- The product should be injected into the dermis or the lower part of the dermis.
- To remove air from the interior of the syringe, press the rod slowly until there is a small drop at the tip of the needle.
- Aspiration before use is recommended when using a needle. Inject content while slowly pulling the needle backwards.
- Stop injection before fully pulling out the needle.
- Use a new needle and cannula for each injection site. Massage the injection site following injection to facilitate conformation to surrounding tissue.
- Periodic additional treatment helps to achieve desired results. The treatment can be adjusted depending on the user's discretion and the preference of the patient.

3. Storage Following Use

- Dispose of used and left over products including needles and cannulas after procedure.
- Due to possible risks of contamination, distortion, and infection, left over products are not to be reused.
- The reuse of this product, a single use medical device, is prohibited.

Precautions During Use

1. Contraindication

- The product is to be used only for intradermal injections.
- This product is not to be injected into blood vessels. Unintended introduction into blood vessels could cause occlusion, ischemia, tissue necrosis, etc.
- When injected into a blood vessel, this may result in serious side effects such as blindness, it is suggested not to use this filler around the

- middle of forehead having a thin skin and a high probability of being injected into blood vessels. So pay special attention when operating.
- Stop injection if the injection site becomes pale, and massage the treatment site until it returns to a normal color.
- The product is not to be used on patients with bleeding disorders or those taking thrombolytic agent or anticoagulants.
- This product is not to be resterilized.
- This product is not to be mixed with other products before injection.
- This product contains lidocaine and should not be used for patients with lidocaine allergy.
- This product is not intended to be used in patients with known hypersensitivity to lidocaine or anesthetic of the amide type.

2. Precautions Regarding Age, Sex, and Physical Conditions

- Do not inject the product into an area where another manufacturer's product is present.
- Do not inject the product into areas where non injectable implants are already present.
- Following injection, the patient should minimize exposure to excessive sunlight or cold until the initial edema and erythema become less severe.
- There is a risk of inflammation of the injection site should the injection area be subject to laser treatment, chemical peeling, or any other procedure based on dermal response before the treatment site is fully healed. The use of this product on an area treated with the previous procedures before full healing of the site can cause inflammation.
- An excessive or too shallow injection can result in the temporary formation of a lump.
- Safety & Effectiveness for Lip Augmentation are not established.
- Do not use this product on pregnant or breastfeeding women or children.

3. Precautions for Undesirable Effects

- The medical practitioners must inform the patients that potential side effects associated with use of this product can occur immediately or may be delayed.
- These include: Unintended injection into blood vessels could cause blindness, occlusion, ischemia, tissue necrosis, stroke, etc. Frequency of 1/10,000~20,000. Edema, bruising, erythema, light pain, tenderness Frequency of 1/50,000: Infection, inflammation, blanching, papules, nodules Low Frequency of the following: itching, hypersensitivity, reactivation of latent for subclinical herpes infection, acne, granuloma, blistering, clotting, facial edema, rash, dermatitis, scarring, skin atrophy, reduction of duration time, ischemia, tissue damage, and augmentation of the capillary. Cases of visual impairment from unintended injections into facial veins have been reported. Cases of inflammation accompanying red spots, edema, tenderness, clotting, and other symptoms have been reported. Cases of pigmentation after inflammation due to melanin accumulation have been reported.
- Safety for a patient being sensitive to keloid pathogenesis, hyperpigmentation and hypertrophic scar is not established.
- The infection site should be treated with adequate treatment or with the removal of the implant. If the symptom is severe, corticosteroid treatment could be effective.

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심볼 / Symbols

| | | | |
|---------------|---------------------------|--|----------------------------------|
| | Catalogue number | | CE Mark |
| | Do not reuse | | Batch code |
| | Caution | | Date of manufacture |
| | Sterilized using steam | | Use by |
| | Temperature limitation | | Fragile, handle with care |
| | Do not resterilize | | Manufacturer |
| | Keep away from sunlight | | Do not use if package is damaged |
| EC REP | Authorized representative | | |

6. Absorbable information

- MONALISA filler is bio-absorbable which means that they will gradually be absorbed into the body. As the gel breaks down, water takes its place, and when totally absorbed, the MONALISA Lidocaine gel disappears unnoticed from the body.
- How long the treatment holds its effect is very individual, but will typically last six to twelve months. This depends on many factors, such as your age, skin type, lifestyle, and muscle activity, as well as on the injection technique.
- There could be some delayed side effects for absorbable of injectable dermal fillers such as inflammation, erythema or edema.

Mass / Packing Unit

Refer to bottom of the package

Storage

1. Store at room temperature (2 ~ 25 ℃)
2. Protect from sunlight
3. Do not freeze the product

Expiration Date

Refer to bottom of the package