# GENOSS Introducer Sheath

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|  |                                  |  |   |
|--|----------------------------------|--|---|
|  | Manufacturer                     |  | Keep away from sunlight   |
|  | Date of manufacture              |  | Keep dry  |
|  | Use by date                      |  | Temperaturd limit   |
|  | Batch code                       |  | Do not re-use   |
|  | Catalogue number                 |  | Consult instructions for use                                    |
|  | Sterilized using ethylene oxide  |  | Caution   |
|  | Do not use if package is damaged |  | Medical device  |
|  | Unique device identifier         |  | Single Sterile barrier system with protective packaging outside |

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GENOSS Introducer Sheath

## B Description

This product is a medical device mainly used to minimize vascular damage that may occur in the process of percutaneously inserting a catheter into a vein or artery during diagnosis and treatment of vascular diseases, and to reduce the inconvenience and complications of the procedure.

## C Intended Use

It is a device that maintains hemostasis while performing interventional procedures and facilitates percutaneous insertion of various devices into veins or arteries.

It is a composite medical device, which consists of the following components.

1. Sheath : sheath used to insert the catheter through the skin.
2. Dilator : manual device that dilates or measures blood vessels.
3. Guide wire : coil-type device inserted into a catheter to guide the catheter into the blood vessel.
4. IV Catheter : tube-shaped catheter that is inserted into a patient's blood vessel for a short period of time (within 30 days) and is used for blood collection and blood pressure monitoring.
5. Introducer needle(Seldinger needle) : disposable needle used for the purpose of healing by inserting a fine needle into the skin or tissue.
6. Blades : It is a manual medical knife used for cutting and incision of human tissue and is disposable.

## D Warnings

1. Be familiar with the product manual before using it.
2. Do not use the product if the sterile packaging of the product is damaged.
3. This product is a disposable sterilization product, so do not re-sterilize or reuse it.
4. Dispose of expired products without using them.
5. It must be used for the using purpose of the product and should not be used for any other purpose.
6. Use immediately after the package opened.

## E Cautions


1. Check if the product is bent or damaged by appearance before use and do not use damaged products.
2. Do not use if the package is opened or the label is damaged and cannot be identified.
3. This device is intended for use only by physicians skilled in percutaneous and intravascular techniques and procedures.
4. Do not apply physical impact to the distal tip part of the catheter to prevent deformation.
5. As this is the disposable sterile medical device, so cannot be reused or re-sterilized.

## F Potential Complications

The following adverse events may include, but are not limited to:

- Vascular trauma (e.g. perforation, dissection)
- Vascular complications (e.g., puncture site complications)
- Blood clots
- Infection/Pain in puncture site
- Arterial spasm
- Embolism/Obstruction
- Bleeding
- Nausea/Vomiting
- Fever/Chills
- Arteriovenous fistula
- Behavioral disorder
- Kidney failure
- Arrhythmia
- Hematoma
- Allergies

## G Preparations for use

Carefully check the GENOSS Introducer sheath to ensure that it is not damaged during handling.

Ensure that the size of the product to be used is suitable for the access vessel and the interventional device used.

Remove Tyvek wrappers and check for all components (Sheath, Dilator, Guide wire, IV. catheter, Introducer needle, Blade).

## H Instructions for use

1. Check and prepare the insertion area using an appropriate aseptic method and local anesthesia if necessary.
2. Remove the GENOSS Introducer sheath and components from the package using a suitable aseptic method.
3. Flushing all components with saline water and wet the surface of the sheath.
4. After flushing using the three-way valve of the sheath, close the valve to maintain the flushing state.
5. A dilator is inserted through the hemostatic valve of the sheath to fix it so that it clicks.
6. Insert a suitable needle into the blood vessel.
  - a) When using an introducer needle, insert the flexible tip of the wire into the blood vessel through the hub of the introducer needle while holding it.
  - b) When using an I.V. catheter, remove the metal needle tube after puncture and insert the flexible tip of the wire into the vessel through a plastic catheter.
7. Position the guide wire in the vessel and remove the needle. After removing the needle, apply pressure to the puncture area to stop bleeding.
8. Insert the assembled sheath and dilator into the puncture area along the guide wire and place it in the intravascular.
9. Bend the hub of the dilator, separate the assembled part, and remove it from the sheath. (The sheath is located in the blood vessel and carefully removes the guide wire at the same time as the dilator is removed)
10. After removing the dilator, suction is performed from the extension line of the sheath to remove remaining air. After suction, flush through a three-way valve using an appropriate solution. At this point, the valve is kept locked to maintain cleaning on the extension line.
11. Removed
  - a) When clinically indicated, the sheath shall be removed. Intravascular compression on the perforated area begins with the sheath slowly being removed. Dispose of sheath as appropriate.
    - ※ When removing the sheath, apply pressure to the puncture area and press it to block the flow of blood as much as possible and remove it slowly.
  - b) When the sheath is left, it is inserted into the sheath using an appropriate size occlusion device. After washing the sheath, the occlusion device is placed in an appropriate position through the sheath and snap. Clean with an appropriate solution through an extension line.
    - ※ Sheath shall be removed within 24 hours of insertion.
12. Dispose of all components including used sheath without reuse.

## H Storage

Store at room temperature (1°C to 30°C) in a cool, dry and dark place.

## I Shelf life

3 years from manufacturing date.

## J Number of uses

This product is a disposable medical device, so it should be not re-sterilized, not reused, and discarded according to the rules after use.

## K Packaging unit

5ea / 1Box

## L Symbols

|  |                                  |  |   |
|--|----------------------------------|--|---|
|  | Manufacturer                     |  | Keep away from sunlight   |
|  | Date of manufacture              |  | Keep dry  |
|  | Use by date                      |  | Temperature limit   |
|  | Batch code                       |  | Do not re-use   |
|  | Catalogue number                 |  | Consult instructions for use                                    |
|  | Sterilized using ethylene oxide  |  | Caution   |
|  | Do not use if package is damaged |  | Medical device  |
|  | Unique device identifier         |  | Single Sterile barrier system with protective packaging outside |

MS-P2759/IFU-DS1902(Rev.2,2603)