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4

English

3

4

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5

6

Intensity control	한국어
	- SONIQ Toothbrush has three modes: low, medium and high. Use it by setting the intensity that suits you. - To set the intensity, give a very short press to the power button while the product power is on. It repeats and changes from low to high. - The SONIQ Electric Toothbrush has a memory function. It remembers the last intensity you used and starts brushing at that intensity the next time you use it.
Toothbrush cleaning and storage How to clean and store your toothbrush head  Toothbrush cleaning and storage 	- Immediately after use, clean the toothbrush head under running water and dry it. - If there is toothpaste left on the toothbrush or if it is not dried properly, it may cause bacterial growth. - After washing and drying, store in a dry, cool place.
5	

Charging	English
	- It can be used for up to 50 days when fully charged. (By intensity: high: 30 days, medium: 40 days, low: 50 days/ based on brushing your teeth for 2 minutes three times a day) The time may vary depending on user. - This product has wireless charger. Gently place the main body on the cradle connected to the power source to start charging. - The light is red if the battery is less than 80%, blue if the battery level is higher than 80% when charging. While in use, the light turns redy when the battery is less than 20%, blue if the battery is higher than 80%. - The toothbrush doesn't work close to the charging cradle for safety. - The charging standard of this product is 5V 1A. Please use an adapter of the appropriate standard.(Charging may not be possible with an adapter less 5V 1A.) - When connecting the USB power source, keep it dry. It may cause fire or electric shock.
Warranty and support Warranty This product warrants for 6 months. After 6 months, costs will be charged. Please contact the salesperson or manufacturer for warranty. (However, even during the warranty period, defects due to the user's negligence, it will void this warranty.)	
6	

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MS-P2886/1FU-DE1901(Rev.3,2606)

FAQ	English
The vibration strength is weak. → If you brush your teeth by pressing the toothbrush hard, the vibration may be weakened. Do not press the toothbrush too hard, but use it soft touch. → If the toothbrush head is connected too strongly , the internal motor may be damaged . Check that there is a gap of about 1mm between the toothbrush head and body . If not , please contact for customer service center. It doesn't work even when charging. → The battery needs to be fully charged for the first time. And the power ON/OFF may be inexperienced. Please read the manual carefully and use it after charging. If it still doesn't work, ask for customer Service centre. → The toothbrush is may be running out of battery . Please use it after enough charging it. If charging does not work even when charged , request customer service center Charging is not working well. → If there are water or dust on the USB connector, contact faulty may occur. Using the USB connector after cleaning and drying. → The wireless charging cradle and the body may not come into contact. Check if there are any foreign substances in the wireless charging cradle, and if charging does not work even after removing it, request replacement. The vibration is too strong. → The vibration may be felt strongly in sensitive and weak areas on the teeth . Use strength as a low mode, or stop using it and take medical attention. Occur the heat while charging. → Due to the characteristics of wireless charging, some amount of heat may exist. If heat still occurs, immediately turn off the charger and contact customer service. There is a strange noise when operating. → The product that uses vibration, a certain amount of noise exists. If the noise continues to occur, request customer service. Is it okay to keep a fully charged toothbrush standing on the charger? → This product that charge finishes after a full charge because there is including prevents overcharging. For the life of the battery, it is recommended to disconnect it from the charger after finishing charge.	
8	

중전	English
Is it okay to leave the charger and toothbrush in a moist place? →The toothbrush is not damaged because of provides IPX8 waterproof. However, it is not recommended as the charger does not support waterproofing. Store in a cool, dry place as much as possible. Is it okay to wash the toothbrush body under running water? → This product is possible because it supports IPX8 waterproof rating.The charger is not waterproof, so do not wash it with water.	
Product information Product name: SONIQ Model number: BSET-001,BSET-02,BSET-03,BSET-04,BSETC-001 Material: ABS, PC, POM and others Rated input: DC 5V, 1A Rated output(Power consumption): 2W Vibration: 30,000/min Wireless charging frequency: 70 KHz ± 10% Batteries: 3.7V, 1100mAh Waterproof rate: IPX8 Dimensions: 180 X 28.5 X 27.5 mm(body) Weight: 116g	
Manufacturer information Manufacturer: genoss co.,Ltd.(Korea) D-factory 56, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do Tel. +82-31-888-5100 www.genoss.com Customer Service The warranty period is 6 months , and even during the free warranty period, it is converted to a paid service in case of repair due to consumer negligence. Tel. +82-31-888-5100 Fax. +82-31-8019-9123	
9	